


# NYLON


(MONOFILAMENT)



NYLON

*AILEE SUTURE*  
ISO 13485 CE 0123

# **Monofilament Nylon**



*The superior smooth surface ensures the effortless passage through the tissue and it is too slippery for bacteria to grow on.*

*Extruded to provide a strong uniform strand, good handling and pliability characteristics make knot easy and secure.*

*As with other synthetic sutures, knot security requires the standard technique of square ties with additional throws.*

*Nylon sutures retain their tensile strength during the healing period with minimal tissue reaction.*

**Non-Absorbable**

**Range gauge sizes 11-0 to 2 (U.S.P)**






















**Blue, Black and White**

**Sterilized by Ethylene Oxide Gas**



























# NYLON

(Monofilament)

Needle		Suture	U.S.P (E.P)				Surgical Use
Length	Description	Color & Length	11-0 (0.1)	10-0 (0.2)	9-0 (0.3)	8-0 (0.4)	
<b>3m/m</b>							
3/8	   <b>75 Microns</b>	 Black 15cm	NK 1102	NK 1002			Micro
		 Black 15cm					
<b>3.5m/m</b>							
3/8	  <b>75 Microns</b>	 Black 15cm	NK 1104	NK 1004			Micro
<b>3.75m/m</b>							
3/8	  <b>75 Microns</b>	 Black 15cm		NK 1041			Micro
<b>4m/m</b>							
3/8	  <b>75 Microns</b>	 Black 15cm	NK 1107				
3/8	  <b>100 Microns</b>	 Black 15cm		NK 1007	NK 907		Micro
		 Black 15cm		NK 1008	NK 908	NK 808	
3/8	  <b>150 Microns</b>	 Black 15cm		NK 1040DN			
<b>4.5m/m</b>							
3/8	  <b>130 Microns</b>	 Black 30cm				NK 812DN	Micro
<b>4.75m/m</b>							
1/2	  <b>100 Microns</b>	 Black 15cm		NK 1042			Micro
3/8	  <b>130 Microns</b>	 Black 15cm		NK 1031	NK 9151	NK 8151	

























# NYLON (Monofilament)

Needle		Suture	U.S.P (E.P)					Surgical Use
Length	Description	Color & Length	11-0 (0.1)	10-0 (0.2)	9-0 (0.3)	8-0 (0.4)	7-0 (0.5)	
<b>5m/m</b>								
3/8			Black 15cm		NK 1051	NK 951		
			Black 15cm		NK 1054	NK 952		
<b>100 Microns</b>								
3/8			Black 15cm		NK 1052			
<b>75 Microns</b>								
<b>Compound Curve</b>								
			Black 30cm		NK 1053DN			
<b>150 Microns</b>								
3/8			Black 15cm		NK 1055	NK 915	NK 815	Cornea, Cataract
			Black 15cm			NK 916	NK 816	
<b>150 Microns</b>								
3/8			Black 30cm				NK 818	
			Black 30cm				NK 817	
			Blue 30cm				NB 817	
<b>300 Microns</b>								
1/2			Black 15cm		NK 1050			
<b>220 Microns</b>								
<b>6m/m</b>								
3/8			Black 30cm		NK 1013			
			Black 30cm	NK 1190DN	NK 1090DN	NK 990DN		
<b>150 Microns</b>								
3/8			Black 30cm		NK 1070	NK 970		Cornea, Cataract
			Black 30cm		NK 1070DN	NK 970DN		
<b>220 Microns</b>								













# NYLON

(Monofilament)

Needle		Suture	U.S.P (E.P)					Surgical Use
Length	Description	Color & Length	10-0 (0.2)	9-0 (0.3)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	
<b>6m/m</b>								
3/8			Blue 30cm			NB 821	NB 721	
			Black 30cm			NK 821	NK 721	
			Blue 30cm			NB 822	NB 722	
			Black 30cm			NK 822	NK 722	
			Black 30cm			NK 821DN		
<b>300 Microns</b>								
1/2			Black 30cm	NK 1080DN				Cataract, Cornea, Vascular
<b>150 Microns</b>								
1/2			Black 30cm	NK 1060				
			Black 30cm	NK 1091DN				
<b>220 Microns</b>								
1/2			Blue 30cm			NB 819	NB 719	
			Black 30cm			NK 819		
<b>300 Microns</b>								
<b>6.5m/m</b>								
3/8			Black 15cm	NK 1065		NK 890		Micro
			Black 30cm	NK 1071DN				
<b>130 Microns</b>								
<b>7m/m</b>								
3/8			Blue 50cm			NB 823	NB 723	
			Black 50cm			NK 823DN		
			Blue 50cm			NB 824	NB 724	Cornea
			Black 50cm			NK 824	NK 724	
			White 50cm			NW 824	NW 724	

# NYLON





























(Monofilament)

Needle		Suture	U.S.P (E.P)					Surgical Use
Length	Description	Color & Length	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	
<b>8m/m</b>								
3/8	 <b>220 Microns</b>	⊙ Black 50cm	NK 825PDN					
		⊙ Blue 50cm	NB 825	NB 725	NB 607			
3/8	 	⊙ Blue 50cm	NB 825DN		NB 607DN			
		⊙ Black 50cm			NK 607			
		⊙ Black 50cm	NK 825DN					
		▽ Blue 50cm	NB 826	NB 726	NB 608			
		▽ Black 50cm	NK 826	NK 726				
		▽ White 50cm	NW 826	NW 726	NW 608			
		▽ Black 50cm	NK 826DN					
3/8	 	▽ Blue 50cm		NB 726P				
		▽ Black 50cm		NK 726P				
		▽ White 50cm		NW 726P				
1/2	 <b>220 Microns</b>	▽ Black 50cm	NK 826P					
3/8	 	▽ White 50cm			NW 608P	NW 508P		
		▽ Blue 50cm			NB 608P	NB 508P		
		▽ Blue 50cm				NB 508PDN		
<b>9m/m</b>								
3/8	 	⊙ Blue 50cm		NB 727	NB 609			
		▽ Blue 50cm		NB 728	NB 610			
		▽ Blue 50cm		NB 728DN	NB 610DN			
		▽ Black 50cm		NK 728	NK 610			
		▽ White 50cm		NW 728	NW 610			
								Plastic, Skin, Vascular




























# NYLON

(Monofilament)

Needle		Suture	U.S.P (E.P)					Surgical Use	
Length	Description	Color & Length	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)		
<b>9m/m</b>									
3/8			Black 50cm			NK 6101P	NK 5101P		Plastic, Skin, Vascular
1/2			Blue 30cm	NB 829	NB 729	NB 611			
<b>11m/m</b>									
3/8			Black 50cm		NK 717P	NK 617P	NK 517P	NK 417P	Plastic, Skin, Vascular
			Blue 50cm		NB 717P	NB 617P	NB 517P		
			White 50cm		NW 717P	NW 617P			
3/8			Blue 50cm		NB 717	NB 617	NB 517		
			Black 50cm		NK 717	NK 617	NK 517		
			White 50cm		NW 717	NW 617	NW 517		
			Blue 50cm			NB 616	NB 516		
			White 50cm			NW 616	NW 516		
			White 50cm				NW 516DN		
			Blue 50cm		NB 717DN	NB 617DN			
3/8			Blue 50cm		NB 717R	NB 617R	NB 517R		
			White 50cm		NW 717R				
1/2			Blue 30cm		NB 731	NB 613	NB 513	NB 413	
			Black 50cm		NK 715P				
1/2			White 50cm			NW 615P	NW 515P		
1/2			Blue 50cm			NB 615			
			Blue 50cm				NB 514	NB 414	
			White 50cm				NW 514		

# NYLON
























(Monofilament)

Needle		Suture	U.S.P (E.P)					Surgical Use	
Length	Description	Color & Length	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)		
<b>13m/m</b>									
3/8		 Blue 50cm			NB 621	NB 521	NB 421		
		 Black 50cm			NK 621	NK 521	NK 421		
		 White 50cm			NW 621	NW 521	NW 421		
		 Black 50cm			NK 621DN				
		 White 50cm			NW 621DN				
		 Blue 50cm				NB 520			
3/8		 Blue 50cm			NB 621P	NB 521P	NB 421P	Plastic, Skin,	
		 Black 50cm				NK 521P			
		 White 50cm			NW 621P	NW 521P	NW 421P		
1/2		 Blue 50cm			NB 619	NB 519			
		 Blue 50cm				NB 518			
1/2		 White 50cm			NW 619 P	NW 519 P			
<b>14m/m</b>									
1/2		 Blue 30cm			NB 622	NB 522	NB 422	NB 300	Vascular
<b>15m/m</b>									
3/8		 Blue 50cm			NB 626	NB 526	NB 426	NB 304	Plastic, Skin
		 Black 35cm			NK 626		NK 426		
		 Black 50cm				NK 526			
		 White 50cm				NW 526	NW 426		
		 Blue 50cm				NB 525	NB 425		



# NYLON

(Monofilament)

Needle		Suture	U.S.P (E.P)					Surgical Use	
Length	Description	Color & Length	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)		
<b>15m/m</b>									
3/8			Black 50cm	NK 626P	NK 526P				Plastic, Skin
			Blue 50cm		NB 526P				
1/2			Blue 50cm		NB 524	NB 424			
<b>17m/m</b>									
3/8			Black 50cm	NK 628P					Plastic, Skin
			Blue 50cm		NB 528P				
3/8			Blue 50cm		NB 528	NB 428	NB 306	NB 206	Plastic, Skin
1/2			Blue 50cm		NB 527	NB 427	NB 305	NB 205	
			Blue 50cm		NB 529	NB 429	NB 307	NB 207	
<b>18m/m</b>									
3/8			Blue 50cm	NB 634	NB 534	NB 434	NB 312	Plastic, Skin, Gas-Int	
			Black 50cm	NK 634	NK 534	NK 434	NK 312		
			White 50cm			NW 434	NW 312		
			White 50cm	NW 634DN					
3/8			Blue 50cm			NB 434P	NB 312P		
1/2			Blue 50cm	NB 630	NB 530	NB 430	NB 308		
			Blue 50cm				NB 310		

# NYLON






(Monofilament)

Needle		Suture	U.S.P (E.P)					Surgical Use	
Length	Description	Color & Length	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)		2-0 (3)
<b>19m/m</b>									
1/2			Blue 50cm		NB 535	NB 435	NB 313		Vascular
1/2			White 50cm	NW 735P					Plastic
			White 50cm	NW 735PDN					
<b>21m/m</b>									
3/8			Blue 50cm		NB 541	NB 441	NB 319	NB 219	
3/8			Blue 50cm			NB 439	NB 317	NB 217	Plastic, Skin, Dental
1/2			Blue 50cm		NB 540	NB 440	NB 318	NB 218	
			Blue 50cm		NB 542	NB 442	NB 320	NB 220	
<b>23m/m</b>									
1/2			Blue 75cm		NB 543	NB 443	NB 321	NB 221	Gas-Int
<b>24m/m</b>									
3/8			Blue 45cm	NB 7024DN					Plastic
			Black 45cm	NK 7024DN					
			White 45cm	NW 7024DN					





















# NYLON

(Monofilament)

Needle		Suture	U.S.P (E.P)					Surgical Use	
Length	Description	Color & Length	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)		
<b>24m/m</b>									
1/2		○	Blue 50cm		NB 544	NB 444	NB 322		Vascular
<b>25m/m</b>									
3/8		▽	Black 50cm	NK 646					
		▽	Blue 50cm		NB 546	NB 446			
		▽	Blue 75cm				NB 324	NB 224	
		▽	Black 75cm				NK 324		
3/8		▽	Blue 50cm		NB 546P		NB 324P		Plastic, Skin, Dental
		▽	Blue 75cm				NB 3241P	NB 224P	
1/2		▽	Blue 50cm					NB 223	
		○	Blue 75cm			NB 447	NB 325	NB 225	
<b>28m/m</b>									
1/2		○	Blue 75cm					NB 226	Gas-Int



















# NYLON

(Monofilament)

Needle		Suture	U.S.P (E.P)					Surgical Use	
Length	Description	Color & Length	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)		
<b>29m/m</b>									
3/8			Blue 75cm	NB 328	NB 228	NB 110	NB 0110	NB 0210	Fascia, Skin
1/2			Blue 75cm		NB 227	NB 109			
			Blue 75cm		NB 229	NB 111	NB 0111		
			Blue 100cm		NB 2291				
1/2 Loop			Blue 150cm		NB 229L				
1/2 Loop			Blue 150cm			NB 111L			
<b>31m/m</b>									
1/2			Blue 75cm		NB 231				Ob-Gyn
<b>33m/m</b>									
3/8			Blue 75cm	NB 333	NB 233	NB 115	NB 0115		Skin, Fascia
1/2			Blue 75cm		NB 234				
1/2 Loop			Blue 150cm				NB 0116L		


















# NYLON (Monofilament)

Needle		Suture	U.S.P (E.P)					Surgical Use
Length	Description	Color & Length	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)	
<b>37m/m</b>								
3/8		 Blue 75cm	NB 339	NB 239	NB 121	NB 0121		Skin, Fascia
		 Black 75cm		NK 239				
1/2		 Blue 75cm		NB 238		NB 0120		
1/2		 Blue 75cm			NB 122	NB 0122		
1/2 Loop		 Blue 150cm				NB 0122L		
<b>41m/m</b>								
3/8		 Blue 75cm		NB 243	NB 125	NB 0125	NB 0225	Skin, Fascia Ob-Gyn
1/2		 Blue 75cm			NB 126	NB 0126		
		 Blue 75cm			NB 124	NB 0124		
1/2		 Blue 75cm			NB 128	NB 0128		
1/2 Loop		 Blue 150cm				NB 126L		

# NYLON












(Monofilament)

Needle		Suture	U.S.P (E.P)					Surgical Use
Length	Description	Color & Length	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)	
<b>45m/m</b>								
3/8		 Blue 100cm		NB 249				
		 Blue 75cm				NB 131	NB 0131	NB 0231
1/2		 Blue 75cm					NB 0132	Skin, Fascia
<b>46m/m</b>								
1/2		 Blue 75cm						NB 0233 Fascia
<b>49m/m</b>								
1/2		 Blue 100cm					NB 0150R	
		 Blue 100cm					NB 0150	
1/2 Loop		 Blue 150cm					NB 150L	Fascia
1/2 Loop		 Blue 150cm					NB 0150L	
		 Blue 200cm					NB 0151L	



# NYLON

(Monofilament)

Needle		Suture	U.S.P (E.P)					Surgical Use	
Length	Description	Color & Length	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)		
<b>55m/m</b>									
			Blue 100cm		NB 255	NB 155			Skin, Subcuticular
<b>60m/m</b>									
			Blue 75cm	NB 360	NB 260				Skin, Subcuticular
	850 Microns								
$\frac{3}{8}$			Blue 100cm		NB 261				Skin, Retention
<b>93m/m</b>									
$\frac{3}{8}$			Blue 100cm		NB 293	NB 193	NB 0193	NB 0293	Skin, Retention
			Blue 75cm				NB 0193DN		
$\frac{3}{8}$	<b>MULTI</b> each with 40mm silicone tubing		Blue 50cm				NB 0193M		
<b>2 Sutures per pack</b>									



## **AILEE COMPANY LIMITED**

Head Office : 280-2 Hakjang-Dong, Sasang-Ku, Pusan, Korea.  
Tel. (051) 322-6161(5 Lines) Fax. (051) 323-1543

Seoul Branch : Yunil Bldg 4F, 1443-15, Seocho-Dong,  
Seocho-Ku, Seoul, Korea.  
Tel. (02) 597-1333~4, Fax. (02) 597-1335

Website : <http://www.ailee.co.kr>