



SILK
(SILICONE TREATED BRAIDED)



AILEE SUTURE
ISO 13485 CE 0123

Braided Silk



A specially developed degumming process removes extraneous material from the raw silk.

This process ensures a compact braid while allowing the filaments to retain their natural body and elasticity.

The braiding on the twisted core makes tensile strength higher. Special treatment for the suppleness and the uniform surface ensures less trauma in the tissue by easy passage-through. Smooth surface helps knot security.

Braided Silk does not become limp or brittle.

The multi-strand presentations may be of benefit where an interrupted technique requires large numbers of sutures.

Non-Absorbable



































Range gauge sizes 10-0 to 3 (U.S.P)

Black, White and Blue

Sterilized by Ethylene Oxide Gas






























SILK

(Silicone Treated Braided)

Needle		Suture	U.S.P (E.P)					Surgical Use
Length	Description	Color & Length	10-0 (0.2)	9-0 (0.3)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	
5m/m								
3/8			Black 30cm			SK 818		Ophthalmic
			Black 50cm				SK 718DN	
	150 Microns							
6m/m								
3/8	Twisted Silk		Blue 30cm			SB 8220		Ophthalmic
			Blue 30cm	SB 1022DN	SB 9222DN	SB 8222DN		
	150 Microns							
3/8			Black 30cm				SK 721	
	220 Microns							
3/8	Twisted Silk		Black 50cm			SK 8221DN		Ophthalmic
			Blue 50cm		SB 9221DN	SB 8221DN		
	300 Microns							
3/8			Black 30cm			SK 822		Ophthalmic
			Black 50cm				SK 722DN	
			Black 50cm				SK 722	
	300 Microns							
3/8			Black 50cm				SK 721	
7m/m								
3/8			Black 30cm			SK 824		Ophthalmic
			Black 50cm				SK 724DN	
			Black 50cm				SK 724	
			Black 30cm			SK 823		
			Black 50cm				SK 723DN	
			Black 50cm				SK 723	
1/2			Black 30cm			SK 8223DN		
	350 Microns							


































SILK

(Silicone Treated Braided)

Needle		Suture	U.S.P. (E.P)					Surgical Use
Length	Description	Color & Length	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	
8m/m								
1/4	 350 Microns	 Black 50cm			SK 6081DN			Ophthalmic Vascular
3/8	 250 Microns	 Black 50cm			SK 6071			
3/8	 220 Microns	 Black 30cm		SK 7251DN				
		 Black 50cm		SK 7261DN				
1/2 Twisted Silk	  220 Microns	 Blue 50cm	SB 8081					
		 Blue 50cm	SB 8081DN					
		 Blue 30cm	SB 8082DN					
3/8	  350 Microns	 Black 50cm	SK 826	SK 726	SK 608			
		 Black 50cm	SK 825	SK 725	SK 607			
		 Black 75cm	SK 8251DN					
		 Black 50cm		SK 726DN	SK 608DN			
9m/m								
3/8	  450 Microns	 Black 50cm		SK 727	SK 609			Ophthalmic Plastic
		 Black 75cm				SK 5171		
		 Black 50cm		SK 728	SK 610			
		 White 50cm			SW 610			
		 Black 50cm		SK 727DN	SK 609DN			
		 Black 50cm			SK 610DN			
1/2	 450 Microns	 Black 30cm		SK 729	SK 611			
1/2	 450 Microns	 Black 30cm				SK 511	SK 411	
3/8	 450 Microns	 White 50cm				SW 5172DN		


































SILK

(Silicone Treated Braided)

Needle		Suture				U.S.P (E.P)				Surgical Use
Length	Description	Color & Length	7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)				
11m/m										
3/8			Black 50cm	SK 716	SK 616	SK 516				
			Black 50cm	SK 717	SK 617	SK 517				
			White 50cm		SW 617	SW 517				
			Black 50cm		SK 617DN	SK 517DN	SK 417DN			
			Black 50cm		SK 616DN					
3/8			Black 50cm		SK 617P	SK 517P	SK 417P	Plastic, Skin, Ophthalmic Vascular		
			Black 50cm		SK 617PDN	SK 517PDN	SK 417PDN			
			Black 50cm	SK 717RDN	SK 617RDN	SK 517RDN	SK 417RDN			
			Black 75cm		SK 617PR					
1/2			Black 50cm	SK 7161						
			Black 30cm				SK 413			
1/2			Black 30cm	SK 731	SK 613	SK 513				
			Black 50cm		SK 614					
			Black 50cm		SK 615					
			Black 50cm		SK 613DN					
350 Microns										
13m/m										
1/2			Black 50cm		SK 619	SK 519		Plastic, Skin, Ophthalmic		
3/8			Black 50cm		SK 621R	SK 521R	SK 421R			
			Black 50cm		SK 621	SK 521	SK 421			
			Black 50cm		SK 621P	SK 521P	SK 421P			
			Black 50cm			SK 520				
3/8			Black 50cm		SK 621DN	SK 521DN	SK 421DN			
			White 50cm				SW 421DN			
3/8			Black 75cm			SK 5211				
14m/m										
1/2			Black 30cm		SK 622	SK 522	SK 422	Vascular		





























SILK

(Silicone Treated Braided)

Needle		Suture	U.S.P (E.P)					Surgical Use
Length	Description	Color & Length	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	
15m/m								
3/8		 Black 50cm		SK 525	SK 425			Ophthalmic Gas-Int, Urol, Skin, Plastic, Vascular
		 Black 50cm	SK 626	SK 526	SK 426	SK 304		
		 Black 75cm		SK 526R	SK 426R			
3/8		 Black 50cm	SK 626P	SK 526P	SK 426P	SK 304P		
3/8	MULTI 4 Sutures per pack	 1 White 3 Black 45cm			SKW 426M SKW 426PM			
1/2		 Black 50cm			SK 424			
1/2		 Black 50cm			SK 424P			
17m/m								
1/2		 Black 50cm		SK 527		SK 305		Ophthalmic Gas-Int, Neuro, Skin, Plastic
		 Black 50cm		SK 529	SK 429	SK 307		
		 Black 75cm				SK 3071		
3/8	 	 Black 50cm		SK 528	SK 428	SK 306		
		 Black 75cm				SK 206		
		 White 50cm			SW 428	SW 306		
3/8		 Black 40cm		SK 528DN	SK 428DN	SK 306DN		
		 Black 50cm		SK 528P	SK 428P	SK 306P		
3/8		 Black 50cm		SK 5291	SK 4291	SK 306R		
1/2	MULTI 6 Sutures per pack	 Black 35cm				SK 307M		
18m/m								
3/8		 Black 50cm		SK 534	SK 434	SK 312		Plastic, Skin, Vascular
3/8		 Black 50cm			SK 434P	SK 312P		
1/2		 Black 50cm		SK 532				
		 Black 50cm		SK 530	SK 430	SK 308	SK 213	























SILK

(Silicone Treated Braided)

Needle		Suture	U.S.P (E.P)					Surgical Use	
Length	Description	Color & Length	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	0 (3.5)		
19m/m									
1/2			Black 50cm	SK 535	SK 435	SK 313			Vascular, Gas-Int
3/8			Black 50cm	SK 5351	SK 4351	SK 3131			
21m/m									
3/8			Black 75cm	SK 5201		SK 3201			Gas-Int, Skin, Dental, Neuro
3/8			Black 50cm		SK 439	SK 317	SK 217		
3/8			Black 50cm	SK 539P	SK 439P	SK 317P			
3/8			Black 50cm		SK 441	SK 319	SK 219		
1/2			Black 75cm	SK 542	SK 442	SK 320	SK 220		
			Black 50cm			SK 318	SK 218		
			Black 55cm		SK 442DN				
1/2	MULTI 		Black 35cm			SK 318M			
	6 Sutures per pack								
1/2			Black 50cm			SK 318P			
23m/m									
1/2			Black 75cm		SK 443	SK 321			Gas-Int
24m/m									
1/2			Black 50cm	SK 544		SK 322			Dental, Gas-Int
25m/m									
3/8	MULTI 		Black 45cm			SK 324M			Skin, Neuro
	4 Sutures Per Pack								
1/2	MULTI 		Black 35cm			SK 325M	SK 225M		
	6 Sutures Per Pack								


SILK

(Silicone Treated Braided)

Needle		Suture	U.S.P. (E.P)					Surgical Use
Length	Description	Color & Length	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	0 (3.5)	
25m/m								
3/8			Black 50cm	SK 546P	SK 446P	SK 324P	SK 224P	
3/8		 	Black 50cm White 75cm	SK 546	SK 446	SK 324	SK 224	
1/2	Removable Needles MULTI	 	Black 45cm Black 75cm		SK 447RM	SK 325RM	SK 225RM	Skin, Oral, Gas-Int, Neuro, Plastic
						SK 3251RM	SK 2251RM	
	8 Sutures per pack							
1/2		   	Black 75cm Black 75cm Black 75cm Black 50cm		SK 447P SK 447	SK 325P SK 325	SK 225P SK 225 SK 223	SK 107P SK 107
5/8			Black 75cm			SK 3241		
28m/m								
1/2			Black 75cm		SK 448	SK 326		Gas-Int
29m/m								
1/2		 	Black 75cm Black 75cm			SK 329 SK 327	SK 229 SK 109	SK 111 SK 109
1/2	MULTI		Black 75cm				SK 229M	Gas-Int Skin
								
	4 Sutures per pack							














SILK

(Silicone Treated Braided)

Needle		Suture	U.S.P (E.P)					Surgical Use
Length	Description	Color & Length	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)	
29m/m								
3/8			Black 75cm	SK 328	SK 228	SK 110	SK0110	Gas-Int, Skin
31m/m								
1/2			Black 75cm			SK112	SK 0112	Ob-Gyn
33m/m								
3/8			Black 75cm	SK 333	SK 233	SK 115	SK 0115	Gas-Int, Skin
1/2		 	Black 75cm Black 75cm	SK 334	SK 234	SK 116	SK 0116 SK 114	
37m/m								
3/8			Black 75cm	SK 339	SK 239	SK 121	SK 0121	Skin
1/2			Black 75cm	SK 340	SK 240		SK 0122	
41m/m								
1/2			Black 75cm		SK 244	SK 126	SK 0126	Gas-Int, Skin, Fascia
3/8			Black 100cm		SK 243	SK 125	SK 0125	
1/2			Black 75cm			SK 0127	SK 0227	
45m/m								
3/8			Black 75cm			SK 131	SK 0131	Skin




SILK

(Silicone Treated Braided)

Needle		Suture	U.S.P (E.P)					Surgical Use
Length	Description	Color & Length	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)	
49m/m								
1/2		 Black 75cm			SK 150			Fascia
55m/m								
		 Black 100cm		SK 255				Skin
MULTI								
		 Black 35cm	SK 355M					Skin
		 Black 45cm	SK 255M					
6 Sutures Per Pack								
60m/m								
		 Black 75cm	SK 360					Skin
	550 Microns							
		 Black 75cm		SK 260	SK 160			Skin
	850 Microns							
3/8		 Black 75cm				SK 0160		

SILK

(Silicone Treated Braided)

Needle		Suture	U.S.P (E.P)					Surgical Use
Length	Description	Color & Length	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)	
75m/m								
		Black 100cm	SK 375	SK 275	SK 175			Skin
80m/m								
$\frac{3}{8}$		Black 75cm	SK 380	SK 280				Skin
90m/m								
$\frac{3}{8}$		Black 100cm		SK 290	SK 190	SK 0190	SK 0290	Skin

Ligature and Non-neededled Sutures

Silk (Silicone Treated Braided Black Silk)

Box of 24 overwrap

USP	1.8m × 1 Lengths	45cm × 10 Lengths	45cm × 17 Lengths	60cm × 13 Lengths	75cm × 10 Lengths
6-0	SK 61801	SK 64510	SK 64517	SK 66013	SK 67510
5-0	SK 51801	SK 54510	SK 54517	SK 56013	SK 57510
4-0	SK 41801	SK 44510	SK 44517	SK 46013	SK 47510
3-0	SK 31801	SK 34510	SK 34517	SK 36013	SK 37510
2-0	SK 21801	SK 24510	SK 24517	SK 26013	SK 27510

Box of 12 overwrap

USP	1.8m × 1 Lengths	45cm × 10 Lengths	45cm × 17 Lengths	60cm × 13 Lengths	75cm × 10 Lengths
0	SK 11801	SK 14510	SK 14517	SK 16013	SK 17510
1	SK 011801	SK 014510	SK 014517	SK 016013	SK 017510
2	SK 021801	SK 024510	SK 024517	SK 026013	SK 027510



AILEE COMPANY LIMITED

Head Office : 280-2 Hakjang-Dong, Sasang-Ku, Pusan, Korea.
Tel. (051) 322-6161(5 Lines) Fax. (051) 323-1543

Seoul Branch : Yunil Bldg 4F, 1443-15, Seocho-Dong,
Seocho-Ku, Seoul, Korea.
Tel. (02) 597-1333~4, Fax. (02) 597-1335

Website : <http://www.ailee.co.kr>