

**AILEE**

SPECIAL STAINLESS STEEL

**SURGICAL  
BLADE**



**AILEE COMPANY LIMITED**

## FOR MINOR OPERATIONS

10



11



12





14




15




# FOR MAJOR OPERATIONS



# SURGICAL BLADES



USE HANDLE NO.3



USE HANDLE NO.4




## Characteristics

1. Sterilized by gamma radiation
2. Made of special stainless steel
3. Grinded keenly with newest auto-grinding machine
4. The best quality guaranteed by special heat-treatment
5. Packed with auto-packing machine
6. International standard goods



# STERILE SURGICAL BLADES



## How to handle goods

For opening aluminium-foil package, grasp both end marked 'open' and pull right and left.



## Sterilized by GAMMA RADIATION

(감마선 소독)  
Sterility guaranteed while package unbroken  
(포장이 파손되지 않는한 소독보장)  
AILEE CO., LTD.  
주식회사 아이리  
280-2, Hakjang-Dong, Buk-ku, Busan, Korea.  
부산시 북구 학장동 280-2



## Aluminium foil packed

Number and shape of blades are imprinted on the box.

For taking out the single peel pack, press the perforated flap at the bottom of 100 blades box.



## Packing unit

1 packet of 100 blades.