

# Coated VICRYL (Polyglactin 910) Suture

## CLASSIFICATION

### Absorbable(흡수성)

가수분해되어지며, 56일에서 70일 사이에 완전히 흡수된다.

### Synthetic(합성사)

90% Glycolide와 10% Lactide의 중합체로 구성된 합성사로 Tissue Reaction이 최소화된다.

### Braided(복선사)

핸들링이 용이하며 매듭의 안전성이 뛰어나다. (10/0나 9/0의 안과용 봉합사 중에는 Monofilament 제품도 있다.)

## COATING

### 코팅재질

Polyglactin 370(Glycolide 30%, Lactide 70%)와 Calcium Sterate로 코팅되어 조직을 부드럽게 통과하고 쉽게 매듭이 내려가도록 재질의 상호작용을 도와준다.

## COLOUR 색상


보라색 또는 흰색


## TENSILE STRENGTH


### 인장강도

상처가 치유되는 기간 동안 강한 힘으로 지탱해주고 그 후에는 빠르게 흡수된다.

% of original strength remaining

적용 14일 후  75%

적용 21일 후  50%

적용 28일 후  25%









적용 35일 후 0%

## INDICATION & CHARACTER

세계적으로 가장 널리 사용되는 봉합사로서, 주로 Soft Tissue의 Approximation이나 안과 수술, 혹은 Ligation 등 대단히 광범위한 부위에 사용된다.

ABSORBABLE  
SYNTHETIC  
BRAIDED

(Polyglactin 910) Suture **Coated VICRYL**

Needle and Description	Length Colour	Size / Product Code						
		10/0	9/0	8/0	7/0	6/0	5/0	4/0
<b>3.9mm</b>								
<b>BV100-3</b> 3/8 Circle 	13cm Violet Monofilament		V402G					
<b>4.2mm</b>								
<b>V130-4</b> TAPER CUT 1/4 Circle 	23cm Violet			J408G				
<b>5.5mm</b>								
<b>CS160-8</b> 1/2 Circle 	30cm Violet Monofilament		V449G					
<b>6.5mm</b>								
<b>CS140-6</b> 3/8 Circle 	10cm Violet Monofilament	V960G						
<b>TG140-8</b> 3/8 Circle 	30cm Violet		V549G	J548G W9560	J566G W9561			
<b>TG100-8</b> 1/4 Circle 	45cm Violet					J544G		
<b>7mm</b>								
<b>CS90-6</b> ULTIMA 1/4 Circle 	10cm Violet Monofilament	V980G						
<b>TG175-8</b> MICRO POINT 1/2 Circle 	45cm Violet				W9565			

COATED VICRYL




















## Coated VICRYL (Polyglactin 910) Suture

ABSORBABLE  
SYNTHETIC  
BRAIDED

Needle and Description	Length Colour	Size / Product Code							
		10/0	9/0	8/0	7/0	6/0	5/0	4/0	
<b>8mm</b>									
<b>S-14</b> 1/4 Circle 		45cm Violet					J570G	J571G	
<b>S-22</b> 1/4 Circle 		45cm Undyed					W9756		
<b>S-24</b> 1/4 Circle 		20cm Violet					W9562		
		30cm Violet					J552G		
		45cm Violet					W9552	W9553	
		45cm Undyed					W9752	W9753	
<b>S-29</b> 1/4 Circle 		45cm Violet					J555G		
<b>G-7</b> <b>PRIME</b> 1/2 Circle 		30cm Violet			W9545				
<b>9.3mm</b>									
<b>BV-1</b> 3/8 Circle 	 Diam. 220 microns	30cm Violet			W9577		W9575		
<b>11mm</b>									
<b>P-1</b> <b>PRIME</b> 3/8 Circle 		45cm Undyed				J488G	W9500T	W9501T	
<b>G-1</b> <b>PRIME</b> 3/8 Circle 		45cm Violet					W9537		
<b>SPAT-10</b> 3/8 Circle 		45cm Undyed					W9760	W9761	W9762

ABSORBABLE  
SYNTHETIC  
BRAIDED

(Polyglactin 910)Suture **Coated VICRYL**

Needle and Description	Length Colour	Size / Product Code						
		9/0	8/0	7/0	6/0	5/0	4/0	3/0
<b>13mm</b>								
<b>TF</b> 1/2 Circle 	 8x45cm Violet						<b>J743D</b> 	
<b>RB-2</b> 1/2 Circle 	 45cm Undyed				<b>W9981</b>	<b>W9982</b>		
<b>C-1</b> 3/8 Circle 	 60cm Violet						<b>W9067</b>	
<b>PC-1</b> <b>PRIME</b> 3/8 Circle 	 45cm Undyed				<b>W9779</b>	<b>W9780</b>		
<b>SC-1</b> <b>STRAIGHT</b> 	 45cm Violet						<b>W9410</b>	
<b>16mm</b>								
<b>CFS-3</b> 3/8 Circle 	 75cm Violet					<b>W9442</b>	<b>W9443</b>	<b>W9444</b>
<b>PC-3</b> <b>PRIME</b> 3/8 Circle 	 45cm Undyed					<b>W9505T</b> <b>W9782</b>	<b>W9506T</b> <b>W9783</b>	<b>W9507T</b> <b>W9784</b>
<b>PS-4C</b> <b>PRIME</b> <b>COMPOUND CURVE</b> 1/2 Circle 	 45cm Undyed					<b>W9436</b>	<b>W9437</b>	










 CONTROL RELEASE

COATED VICRYL



**Coated VICRYL** (Polyglactin 910) Suture

ABSORBABLE  
SYNTHETIC  
BRAIDED

Needle and Description	Length Colour	Size / Product Code						
		7/0	6/0	5/0	4/0	3/0	2/0	1/0
<b>17mm</b>								
<b>BB</b> 3/8 Circle ✓ 	45cm Violet				W9074			
<b>RB-1</b> ✓ 1/2 Circle ✓ 	8x45cm Violet					J713D 		
	8x45cm Undyed				J714D 			
<b>RB-1 PLUS</b> 1/2 Circle 	75cm Violet			W9105	W9106			
<b>UV-17</b> TAPERCUT 5/8 Circle 	75cm Undyed				W9509T			
<b>19mm</b>								
<b>FS-2</b> 3/8 Circle 	45cm Violet				W9386			
	45cm Undyed				W9951			
<b>PC-5</b> PRIME 3/8 Circle 	45cm Undyed				W9510T	W9511T		
	75cm Undyed				W9570T	W9571T		
<b>PS-2</b> PRIME 3/8 Circle 	45cm Undyed				W9515T			
<b>20mm</b>								
<b>SH-2 PLUS</b> 1/2 Circle 	75cm Violet				W9113	W9114		

ABSORBABLE  
SYNTHETIC  
BRAIDED

(Polyglactin 910) Suture **Coated VICRYL**

Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
<b>22mm</b>								
<b>SH-1</b> 1/2 Circle ⊙		8x45cm Violet		<b>J771D</b> <b>CR</b>	<b>J772D</b> <b>CR</b>			
<b>CX-1</b> 1/2 Circle △		45cm Violet			W9472			
<b>ASH-1</b> J NEEDLE ⊙		75cm Violet		W9224	W9223			
<b>SKI-22</b> ⊙		20cm Violet			E9903S <b>VB</b>	E9902S <b>VB</b>		
<b>26mm</b>								
<b>SH</b> 1/2 Circle ⊙		8x45cm Violet			<b>J774D</b> <b>CR</b>			
<b>SH PLUS</b> 1/2 Circle ⊙		75cm Violet			W9120	W9121		
<b>UR-6</b> 5/8 Circle ⊙		70cm Violet				J602H	J603H	
<b>V-7</b> TAPERCUT 1/2 Circle ⊙		75cm Violet			W9180	W9181		
<b>V-7 HEAVY</b> TAPERCUT 1/2 Circle ⊙		75cm Violet				W9350		
<b>PC-25</b> PRIME 3/8 Circle △		75cm Undyed			W9581T	W9582T		
<b>FS</b> 3/8 Circle ▽		75cm Undyed			W9890			

COATED VICRYL

















**CR** CONTROL RELEASE

**VB** VISI-BLACK













**Coated VICRYL** (Polyglactin 910) Suture

 ABSORBABLE  
 SYNTHETIC  
 BRAIDED

Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
<b>31mm</b>								
<b>MO-5</b> 1/2 Circle 		75cm Violet					W9215	
<b>MH-1 PLUS</b> 1/2 Circle 		75cm Violet			W9130	W9136	W9138	W9213
<b>V-30</b> <b>TAPERCUT</b> 1/2 Circle 		75cm Violet				W9360		
<b>36mm</b>								
<b>MH plus</b> 1/2 Circle 		75cm Violet				W9140	W9141	
		90cm Violet				W9440		
<b>V-34</b> <b>TAPERCUT</b> 1/2 Circle 		75cm Violet					W9364	
<b>V-34 heavy</b> <b>TAPERCUT</b> 1/2 Circle 		75cm Violet						W9365
<b>ASH-35</b> <b>J NEEDLE</b> 		75cm Violet				W9220	W9221	
<b>FSLX</b> 3/8 Circle 		75cm Violet				W9390		

ABSORBABLE  
SYNTHETIC  
BRAIDED





(Polyglactin 910) Suture **Coated VICRYL**

Needle and Description	Length Colour	Size / Product Code					
		5/0	4/0	3/0	2/0	1/0	1
<b>40mm</b>							
<b>CT</b> 1/2 Circle  	75cm Violet					W9230	W9231
	90cm Violet					W9430	W9431
<b>LH</b> 1/2 Circle  	75cm Violet				W9150		
<b>V-37</b> <b>TAPERCUT</b> 1/2 Circle  	75cm Violet						W9368
<b>V-37 heavy</b> <b>TAPERCUT</b> 1/2 Circle  	75cm Violet						W9373
<b>CP</b> 1/2 Circle  	75cm Violet					W9320	W9321

COATED VICRYL



Coated VICRYL (Polyglactin 910) Suture

Needle and Description	Length Colour	Size / Product Code				
		3/0	2/0	1/0	1	2
<b>45mm</b>						
<b>V-39 heavy</b> TAPERCUT 1/2 Circle 	75cm Violet				W9377	
<b>48mm</b>						
<b>CTX</b> 1/2 Circle 	75cm Violet				W9251	
	90cm Violet				W9451	
<b>65mm</b>						
<b>TP-1</b> 1/2 Circle 	4x70cm Violet				J650G	
<b>80mm</b>						
<b>MO-80</b> 1/2 Circle 	100cm Violet				W9289	

COATED VICRYL

ABSORBABLE  
SYNTHETIC  
BRAIDED

Coated VICRYL NON - NEEDLED

Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
SUTUPAK* Pre-Cut Sutures for Ligation	12x45cm Violet		V903H	V904H	V905E	V906E	V907E	V908E
	6x45cm Violet	V1223H	V1224H	V1225H	V1226H			
	5x70cm Violet		V1214E	V1215E	V1216E	V1217E	V1218H	V1219H

COATED VICRYL