

Coated

**VICRYL\***  
PlusAntibacterial  
(Polyglactin 910)  
Suture

## COATED VICRYL Plus (Polyglactin 910) Suture

### CLASSIFICATION

#### Absorbable(흡수성)

가수분해 되어지며, 56일에서 70일 사이에 완전히 흡수된다.

#### Synthetic(합성사)

90% Glycolide와 10% Lactide의 중합체로 구성된 합성사로 Tissue Reaction이 최소화된다.

#### Braided(복선사)

핸들링이 용이하며 매듭의 안전성이 뛰어나다.

### COATING

#### 코팅재질

Polyglactin 370(Glycolide 30%, Lactide 70%)와 Calcium Sterate로 코팅되어 조직을 부드럽게 통과하고 쉽게 매듭이 내려가도록 재질의 상호작용을 도와준다.

### COLOUR

#### 색상

보라색 또는 흰색

### TENSILE STRENGTH

#### 인장강도

상처가 치유되는 기간 동안 강한 힘으로 지탱해주고 그 후에는 빠르게 흡수된다.

### INDICATION

Soft Tissue Approximation이나 Ligation 등 대단히 광범위한 부위에 사용한다.

### CHARACTER

Coated Vicryl Plus는 세계 최초의 Antibacterial Suture로서 IRGACARE® MP 라는 초순수 항균인자를 함유하고 있어 Suture 부위에 Bacterial Colonization 현상을 효과적으로 막아준다.

- Coated Vicryl Plus는 미식품의약품안전청(FDA) 승인을 득한 제품이다.
- IRGACARE® MP - A broad spectrum antibacterial agent with proven efficacy
- IRGACARE® MP 는 흔히 창상감염을 유발하는 Staphylococcus Aureus, Staphylococcus Epidermidis, MRSA, MRSE 등 주요 병원균에 효과적으로 대응하는 방어막을 형성하며, 미국에서 30년 이상 Over-the-counter Healthcare Product에 사용 되어 그 안전성이 입증되었다.

% of original strength remaining

적용 14일 후 75%

적용 21일 후 50%

적용 28일 후 25%

적용 35일 후 0%

#### Antibacterial Suture



VICRYL Plus

VICRYL

ABSORBABLE  
SYNTHETIC  
BRAIDED

(Polyglactin 910) Suture **COATED VICRYL Plus**




















COATED VICRYL Plus

Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
<b>8mm</b>								
<b>P-2</b> 1/2 Circle ▽	45cm Undyed	VCP503G	VCP504G					
<b>11mm</b>								
<b>TF-1 PLUS</b> 1/2 Circle ⊙	70cm Violet	VCP133H						
<b>P-1</b> 3/8 Circle ▽	45cm Undyed	VCP490G						
<b>13mm</b>								
<b>RB-2</b> 1/2 Circle ⊙	45cm Undyed	VCP9982H						
<b>C-1</b> 3/8 Circle ⊙	70cm Violet		VCP9067H					
<b>TF</b> 1/2 Circle ⊙	8x45cm Violet		VCP743D <b>CR</b>					
	70cm Undyed	VCP433H	VCP434H					
<b>PC-1</b> <b>PRIME</b> 3/8 Circle △	45cm Undyed	VCP834G	VCP835G					
<b>P-3</b> <b>PRIME</b> 3/8 Circle ▽	45cm Violet	VCP463G	VCP464G					
	45cm Undyed	VCP493G	VCP494G					
<b>PS-5</b> 1/2 Circle ▽	45cm Undyed	VCP593G	VCP594G					
<b>SC-1</b> <b>STRAIGHT</b> ▽	45cm Violet		VCP9410H					



Coated VICRYL Plus (Polyglactin 910) Suture

COATED VICRYL Plus

Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
<b>16mm</b>								
<b>PC-3</b> PRIME 3/8 Circle 	 45cm Undyed	VCP9505H	VCP9506H	VCP9507H				
<b>PS-3</b> 3/8 Circle 	 45cm Undyed	VCP500G						
<b>PS-4</b> 1/2 Circle 	 45cm Undyed	VCP506G	VCP507G					
<b>PS-4C</b> COMPOUND CURVE 1/2 Circle 	 45cm Undyed		VCP656G					
<b>FS-3</b> 3/8 Circle 	 70cm Violet	VCP9442H	VCP394H	VCP395H				
<b>17mm</b>								
<b>RB-1</b> 1/2 Circle 	 70cm Violet	VCP303H	VCP304H	VCP305H	VCP306H			
	8x45cm Violet			VCP713D 				
	8x45cm Undyed		VCP714D 					
	8x70cm Violet		VCP799E 					
<b>BB PLUS</b> 3/8 Circle 	 45cm Violet		VCP9074H					
<b>V-5</b> 1/2 Circle 	 70cm Undyed		VCP230H	VCP2644H				

ABSORBABLE  
SYNTHETIC  
BRAIDED





(Polyglactin 910) Suture **Coated VICRYL Plus**

Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
<b>19mm</b>								
<b>PC-5</b> PRIME 3/8 Circle ▲	45cm Undyed	VCP822G	VCP823G	VCP824G				
	70cm Undyed		VCP9570H	VCP9571H				
<b>PS-2</b> PRIME 3/8 Circle ▼	45cm Undyed	VCP495G	VCP496G	VCP497G				
	70cm Undyed		VCP426H	VCP427H	VCP428H			
<b>22mm</b>								
<b>CT-3</b> 1/2 Circle ⊙	70cm Violet			VCP327H	VCP328H	VCP329H		
<b>MO-7</b> 1/2 Circle ⊙	70cm Violet						VCP250H	VCP251H
<b>SH-1</b> 1/2 Circle ⊙	70cm Violet		VCP310H	VCP311H	VCP312H			
	70cm Undyed		VCP218H	VCP219H				
	8x45cm Violet		VCP771D <b>CR</b>	VCP772D <b>CR</b>	VCP792E <b>CR</b>			
<b>SH-1 PLUS</b> 1/2 Circle ⊙	8x70cm Violet		VCP798G <b>CR</b>	VCP795E <b>CR</b>				
<b>V-6</b> 1/2 Circle ⊙	70cm Violet			VCP277H	VCP278H			



Coated VICRYL Plus (Polyglactin 910) Suture





COATED VICRYL Plus

Needle and Description	Length Colour	Size / Product Code							
		5/0	4/0	3/0	2/0	1/0	1	2	
<b>22mm</b>									
<b>OS-4</b> 1/2 Circle 	3x45cm Violet				VCP706T <b>CR</b>	VCP707T <b>CR</b>	VCP708T <b>CR</b>		
	70cm Violet							VCP9287H	
	70cm Undyed					VCP694H	VCP695H		
<b>X-1</b> 1/2 Circle 	70cm Violet			VCP460H	VCP461H				
	70cm Undyed			VCP458H	VCP459H				
	8x45cm Violet			VCP790D <b>CR</b>					
<b>24mm</b>									
<b>PS-1</b> PRIME 3/8 Circle 	45cm Undyed		VCP682G	VCP683H	VCP684H				
	<b>FS-1</b> 3/8 Circle 	70cm Violet		VCP451H	VCP452H	VCP453H			
70cm Undyed			VCP441H	VCP442H	VCP443H				

ABSORBABLE  
SYNTHETIC  
BRAIDED

(Polyglactin 910) Suture **Coated VICRYL Plus**

COATED VICRYL Plus








Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
<b>26mm</b>								
<b>SH</b> 1/2 Circle 	70cm Violet		VCP315H	VCP316H	VCP317H	VCP318H		
	70cm Undyed		VCP415H	VCP416H	VCP417H	VCP418H		
	90cm Violet		VCP241H	VCP242H	VCP243H	VCP244H		
	8x45cm Violet		VCP773D <b>CR</b>	VCP774D <b>CR</b>	VCP775D <b>CR</b>			
	8x70cm Violet			VCP784D <b>CR</b>	VCP785D <b>CR</b>			
<b>SHB</b> ETHIGUARD BLUNT POINT 1/2 Circle 	8x45cm Undyed		VCPB864D <b>CR</b>					
<b>CT-2</b> 1/2 Circle 	70cm Violet		VCP332H	VCP333H	VCP334H	VCP335H	VCP336H	
	8x45cm Violet				VCP726D <b>CR</b>	VCP727D <b>CR</b>		
	4x45cm Violet						VCP746E <b>CR</b>	VCP756E <b>CR</b>
	8x70cm Violet					VCP747G <b>CR</b>	VCP748G <b>CR</b>	
<b>UR-6</b> 5/8 Circle 	70cm Violet				VCP602H	VCP603H		

**CR** CONTROL RELEASE



Coated VICRYL Plus (Polyglactin 910) Suture

COATED VICRYL Plus








Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
<b>26mm</b>								
<b>MO-6</b> 1/2 Circle 	90cm Violet						VCP237H	VCP238H
	8x45cm Violet				VCP789D <b>CR</b>			
<b>V-7</b> 1/2 Circle 	70cm Violet			VCP998H	VCP999H	VCP279H		
<b>PC-25</b> PRIME 3/8 Circle 	70cm Undyed		VCP9580H	VCP9581H	VCP9582H			
<b>TN</b> TROCAR POINT COMPOUND CURVE 1/2 Circle 	110cm Violet				MIC507H	MIC509H		
<b>CP-2</b> 1/2 Circle 	70cm Undyed			VCP868H	VCP869H	VCP870H	VCP871H	
	8x45cm Undyed			VCP761D <b>CR</b>	VCP762D <b>CR</b>			
<b>30mm</b>								
<b>PSL</b> 3/8 Circle 	70cm Undyed			VCP502H	VCP596H			
<b>FSL</b> 3/8 Circle 	70cm Violet			VCP585H	VCP586H	VCP587H		

**CR** CONTROL RELEASE

ABSORBABLE  
SYNTHETIC  
BRAIDED

(Polyglactin 910) Suture **Coated VICRYL Plus**



COATED VICRYL Plus

Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
<b>31mm</b>								
<b>MH-1 PLUS</b> 1/2 Circle 	70cm Violet			VCP245H	VCP246H	VCP247H	VCP9213H	
	70cm Violet				VCP1036H	VCP1037H	VCP1038H	
	70cm Violet					VCP9215H		
	90cm Violet						VCP253H	
	8x45cm Violet					VCP801D 		
<b>MO-5</b> 1/2 Circle 	8x70cm Violet				VCPP81D 	VCPP80D 		
	70cm Violet				VCP9360H	VCP9361H	VCP9362H	
<b>V-30</b> 1/2 Circle 	70cm Violet						VCP9335H	
<b>V-38</b> 1/2 Circle 	70cm Violet							



## Coated VICRYL Plus (Polyglactin 910) Suture

COATED VICRYL Plus





Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
<b>36mm</b>								
<b>CT-1</b> 1/2 Circle 	70cm Violet			VCP338H	VCP339H	VCP340H	VCP341H	
	90cm Violet			VCP344H	VCP345H	VCP346H	VCP347H	VCP348H
	8x45cm Violet			VCP738D <b>CR</b>	VCP739D <b>CR</b>	VCP740D <b>CR</b>	VCP741D <b>CR</b>	VCP742G <b>CR</b>
	8x70cm Violet					VCPP31D <b>CR</b>		
	4x70cm Violet						VCP736G <b>CR</b>	VCP737G <b>CR</b>
	8x70cm Undyed				VCPP42D <b>CR</b>	VCPP41D <b>CR</b>	VCPP40D <b>CR</b>	
<b>CTB-1</b> ETHIGUARD BLUNT POINT 1/2 Circle 	70cm Violet					VCPB340H	VCPB341H	
	70cm Undyed				VCPB259H			
	90cm Undyed			VCPB944H	VCPB945H	VCPB946H	VCPB947H	
	8x45cm Undyed				VCPB839D <b>CR</b>	VCPB840D <b>CR</b>	VCPB841D <b>CR</b>	

**CR** CONTROL RELEASE

ABSORBABLE  
SYNTHETIC  
BRAIDED

(Polyglactin 910) Suture **Coated VICRYL Plus**

COATED VICRYL Plus






Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
<b>36mm</b>								
<b>MH PLUS</b> 1/2 Circle 	70cm Violet			VCP322H	VCP323H	VCP324H	VCP325H	
	90cm Violet			VCP100H	VCP101H	VCP102H	VCP103H	
<b>MO-4</b> 1/2 Circle 	70cm Violet					VCP436H	VCP437H	
	90cm Violet							VCP438H
	8x45cm Violet				VCP700D <b>CR</b>	VCP701D <b>CR</b>	VCP702D <b>CR</b>	
	8x70cm Violet					VCPP71D <b>CR</b>		
<b>UR-5</b> 5/8 Circle 	70cm Violet				VCP375H	VCP376H		
<b>V-34</b> 1/2 Circle 	90cm Violet			VCP516H	VCP517H	VCP518H	VCP519H	VCP520H



Coated VICRYL Plus (Polyglactin 910) Suture




ABSORBABLE  
SYNTHETIC  
BRAIDED

COATED VICRYL Plus

Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
<b>36mm</b>								
<b>ASH-35</b> J NEEDLE 	70cm Violet					VCP9221H		
	70cm Violet				VCP466H	VCP467H	VCP468H	
<b>CP-1</b> 1/2 Circle 	70cm Undyed			VCP265H	VCP266H	VCP267H	VCP268H	
	90cm Violet				VCP472H	VCP473H	VCP474H	VCP482H
	90cm Undyed				VCP1472H	VCP1473H	VCP1474H	VCP1475H
	70cm Undyed				VCP533H	VCP534H	VCP535H	
<b>OS-6</b> 1/2 Circle 	70cm Undyed				VCP533H	VCP534H	VCP535H	
	3x45cm Violet				VCP710T 	VCP711T 		

ABSORBABLE  
SYNTHETIC  
BRAIDED






(Polyglactin 910) Suture **Coated VICRYL Plus**

Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
<b>40mm</b>								
<b>CT</b> 1/2 Circle 	70cm Violet				VCP351H	VCP352H	VCP353H	
	90cm Violet			VCP356H	VCP357H	VCP358H	VCP359H	VCP360H
	150cm Loop Violet						VCP1408E	
	8x45cm Violet				VCP751D <b>CR</b>	VCP752D <b>CR</b>	VCP753D <b>CR</b>	VCP754G <b>CR</b>
	4x70cm Violet						VCP766G <b>CR</b>	VCP767G <b>CR</b>
<b>CTB</b> ETHIGUARD BLUNT POINT 1/2 Circle 	90cm Violet					VCP9996H	VCP9997H	
<b>V-37</b> 1/2 Circle 	70cm Violet						VCP9373H	
	70cm Undyed				VCP9900H	VCP9901H		






**Coated VICRYL Plus** (Polyglactin 910) Suture

**COATED VICRYL Plus**

Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
<b>40mm</b>								
<b>CP</b> 1/2 Circle ▽	 70cm Violet				VCP478H	VCP479H	VCP480H	
						VCP485H	VCP486H	VCP481H
<b>OS-8</b> 1/2 Circle ▽	 70cm Undyed					VCP698H	VCP699H	
					VCP716T <b>CR</b>	VCP717T <b>CR</b>	VCP718T <b>CR</b>	
<b>45mm</b>								
<b>MO-45</b> 1/2 Circle ⊙	 70cm Violet						VCP9245H	VCP9246H
<b>BT-3</b> ETHIGUARD BLUNT POINT 1/2 Circle ⊙	 90cm Violet						VCP9998H	VCP9999H
<b>V-39</b> 1/2 Circle ▽	 70cm Violet				VCP9375H	VCP9376H	VCP9377H	VCP9378H

ABSORBABLE  
SYNTHETIC  
BRAIDED

(Polyglactin 910) Suture **Coated VICRYL Plus**







Needle and Description	Length Colour	Size / Product Code							
		5/0	4/0	3/0	2/0	1/0	1	2	
<b>48mm</b>									
<b>CTX</b> 1/2 Circle ⊙		70cm Violet				VCP363H	VCP364H	VCP365H	
		90cm Violet				VCP369H	VCP370H	VCP371H	VCP372H
		150cm Loop Violet						VCP1428E	
		8x45cm Violet					VCP764D <b>CR</b>	VCP765D <b>CR</b>	VCP755G <b>CR</b>
		8x70cm Violet							VCP2577G <b>CR</b>
<b>CTXB</b> ETHIGUARD BLUNT POINT 1/2 Circle ⊙		90cm Undyed					VCPB978H	VCPB977H	
		8x45cm Undyed					VCPB724D <b>CR</b>	VCPB725D <b>CR</b>	
<b>CPX</b> 1/2 Circle ▽		90cm Violet						VCP1058H	VCP1059H

**CR** CONTROL RELEASE



ABSORBABLE  
SYNTHETIC  
BRAIDED

(Polyglactin 910) Suture **Coated VICRYL Plus**

Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
<b>48mm</b>								
<b>CTX</b> 1/2 Circle  	70cm Violet				VCP363H	VCP364H	VCP365H	
	90cm Violet				VCP369H	VCP370H	VCP371H	VCP372H
	150cm Loop Violet						VCP1428E	
	8x45cm Violet					VCP764D <b>CR</b>	VCP765D <b>CR</b>	VCP755G <b>CR</b>
	8x70cm Violet							VCP2577G <b>CR</b>
<b>CTXB</b> ETHIGUARD BLUNT POINT 1/2 Circle  	90cm Undyed					VCPB978H	VCPB977H	
	8x45cm Undyed					VCPB724D <b>CR</b>	VCPB725D <b>CR</b>	
<b>CPX</b> 1/2 Circle  	90cm Violet						VCP1058H	VCP1059H






**CR** CONTROL RELEASE

COATED VICRYL Plus

**Coated VICRYL Plus** (Polyglactin 910) Suture

ABSORBABLE  
SYNTHETIC  
BRAIDED

COATED VICRYL Plus

Needle and Description	Length Colour	Size / Product Code					
		4/0	3/0	2/0	1/0	1	2
<b>60mm</b>							
<b>KS</b> STRAIGHT 	70cm Undyed	VCP662H	VCP663H	VCP664H			
<b>65mm</b>							
<b>BP-5</b> BLUNT POINT 3/8 Circle 	100cm Violet					VCP9391G	
<b>70mm</b>							
<b>XLH</b> 1/2 Circle 	70cm Violet			VCP581G		VCP583G	
<b>75mm</b>							
<b>CTXX</b> 1/2 Circle 	90cm Violet						VCP368G
<b>80mm</b>							
<b>MO-80</b> 1/2 Circle 	100cm Violet					VCP9289G	



ABSORBABLE  
SYNTHETIC  
BRAIDED

## Coated VICRYL Plus NON-NEEDED

Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	1/0	1	2
SUTUPAK* Pre-Cut Sutures for Ligation	12x45cm Violet		VCP103G	VCP104G	VCP105G	VCP106G	VCP107G	
	6x45cm Violet		VCP1224H	VCP1225H	VCP1226H	VCP1227E	VCP1228E	
	5x70cm Violet		VCP1214E	VCP1215E	VCP1216E	VCP1217E	VCP1218H	