

MERSILK (Silk) Suture

CLASSIFICATION

Non-Absorbable(비흡수성)

비흡수성 봉합사로 약 3개월 이상 오랜 기간 상처를 지탱할 수 있다.

Natural(자연사)

Degumming 공법을 통하여 99% 이상의 Protein으로 구성되어 있는 순도가 높은 천연 Silk로 만들어져 있다.

Braided(복선사)

봉합사 중에서 가장 뛰어난 핸들링과 매듭의 안정성을 지니고 있다.

COATING

코팅재질

Bees Wax로 코팅처리하여 강도를 높이고 복선사의 단점인 모세관 현상을 줄여 Bacteria에 감염될 위험을 최소화 시켰다.

COLOUR

색상








검정색 또는 아이보리색

INDICATION & CHARACTER

일반적으로 피부봉합이나 Ligation에 많이 이용되지만 그 외에도 다양한 용도에 사용된다.

NON-ABSORBABLE
NATURAL
BRAIDED


















(Silk)Suture **MERSILK**

Needle and Description	Length Colour	Size / Product Code						
		8/0	7/0	6/0	5/0	4/0	3/0	2/0
6.5mm								
TG140-8 MICROPOINT 3/8 circle 	30cm Black		W1814					
8mm								
G-6 3/8 circle 	45cm Black		W817					
S-24 1/4 circle 	45cm Black			W888				
9.3mm								
BV-1 3/8 circle 	45cm Black			W593				
11mm								
BV 3/8 circle 	75cm Black				W595			
P-1 PRIME 3/8 circle 	45cm Black			W812				
P-1 PRIME 3/8 circle 	45cm Black			W815				

Needle and Description		Length	Size / Product Code							
			8/0	7/0	6/0	5/0	4/0	3/0	2/0	
13mm										
TF	1/2 circle	8x45cm								CO84D
C-2	CONVENTIONAL CUTTING	75cm								W468H
C-2	3/8 circle	75cm								W580
C-3	3/8 circle	45cm								W811
16mm										
PC-3	PRIME	4x45cm	3 Black	1 Ivory						W606
PC-3	3/8 circle	45cm	Black							W529H
										W528H
										W604
										W605
CFS-3	3/8 circle	4x45cm	3 Black	1 Ivory						W505
										W501H
										W502H
										W503H
CFS-3	3/8 circle	38cm	Black							W641

NON-ABSORBABLE
NATURAL
BRAIDED

(Silk)Suture **MERSILK**







Needle and Description	Length Colour	Size / Product Code						
		8/0	7/0	6/0	5/0	4/0	3/0	2/0
17mm								
BB 3/8 circle 	 75cm Black				W581	W582		
RB-1 3/8 circle 	8x45cm Black					C054D 	C053D 	
	75cm Black					W549H	W546H	
19mm								
FS-2 3/8 circle 	 45cm Black					W329H		
20mm								
BB-2 3/8 circle 	 75cm Black						W598H	
SH-2 1/2 circle 	 75cm Black					W586H	W587H	
22mm								
SH-1 1/2 circle 	 8x45cm Black					C004D 	C003D 	
CX-1 3/8 circle 	 45cm Black						W577H	W578H

MERSILK

Needle and Description		Length	Colour	Size / Product Code					
		4/0	3/0	2/0	0	1	2	3	
26mm									
SH 1/2 Circle		75cm Black	W570H	W327H					
		8x45cm Black	C014D	C013D	C012D				
		8x75cm Black	C018D	C017D	C016D				
		8x45cm Black			C027D				
CT-2 1/2 Circle		8x45cm Black							
		45cm Black	W328H	W321H					
		45cm Black	W537H						
PC-25 PRIME 3/8 Circle		45cm Black							
		45cm Black							
		45cm Black							
31mm									
MH-1 1/2 Circle		75cm Black	W571H	W333H	W334H				
		75cm Black						W723H	
MO-5 1/2 Circle		75cm Black							

NON-ABSORBABLE
NATURAL
BRAIDED

(Silk)Suture **MERSILK**

Needle and Description	Length Colour	Size / Product Code						
		4/0	3/0	2/0	0	1	2	3
36mm								
MH 1/2 Circle 	75cm Black		W583			W338H		
CT-1 1/2 Circle 	8x45cm Black				C021D 			
FSLX 3/8 Circle 	75cm Black			W667H				
40mm								
CT 1/2 Circle 	75cm Black			W584				
45mm								
LS-1 3/8 Circle 	75cm Black				W782H			

MERSILK

