

ETHILON (Polyamid 6 and 66) Suture

CLASSIFICATION

Non-Absorbable(비흡수성)

조직내에서 영구적인 인장강도를 유지한다.

Synthetic(합성사)

Polyamid(Surgical grade of nylon)으로 만들어져 Tissue Reaction이 최소화된다.

Monofilament(단선사)

표면이 매끄러워 조직 투과력이 뛰어나고 박테리아의 증식을 최소화 할 수 있어 감염의 위험이 적다.

COLOUR

색상

















파란색, 검정색 혹은 무색

INDICATION & CHARACTER

단선사이지만 매우 유연하고 핸들링이 뛰어나 주로 피부 봉합에 많이 사용되며, 재질의 특성상 매우 가는 봉합사로도 생산이 가능하여 안과 수술이나 미세 수술에도 많이 사용된다.

NON-ABSORBABLE
SYNTHETIC
MONOFILAMENT

(Polyamid 6 and 66) Suture **ETHILON**



















Needle and Description	Length Colour	Size / Product Code							
		11/0	10/0	9/0	8/0	7/0	6/0	5/0	
3.5mm									
V100-3 MICRO POINT TAPERCUT 3/8 Circle 	 Diam. 100 microns chord length 3.0 mm	13cm Black		W2899	W2898				
3.8mm									
BV50-3 3/8 Circle 	 Diam. 50 microns chord length 3.2 mm	13cm Black	W2881						
BV75-3 3/8 Circle 	 Diam. 75 microns chord length 3.0 mm	13cm Black		W2870	W2871				
4.5mm									
TG135-6 3/8 Circle 		13cm Black				W2901			
V150-4 MICRO POINT TAPERCUT 1/4 Circle 	 Diam. 150 microns	15cm Black		W1701					
4.7mm									
BVH100-3 1/2 Circle 	 Diam. 100 microns	13cm Black		W2860					
5mm									
BV130-4 3/8 Circle 	 Diam. 130 microns chord length 3.7 mm	13cm Black		W2814	W2813 W2913 VB	W2812			
TG6-C ADVANCED MICRO POINT COMPOUND CURVE 	 Diam. 150 microns	30cm Black		W7799					

VB VISI-BLACK

ETHILON

ETHILON (Polyamid 6 and 66)Suture

NON-ABSORBABLE
SYNTHETIC
MONOFILAMENT

Needle and Description	Length Colour	Size / Product Code					
		11/0	10/0	9/0	8/0	7/0	6/0
5.1mm							
BV100-4 3/8 Circle   Diam. 100 microns chord length 4.0 mm	13cm Black		W2830	W2829			
5.5mm							
CS160-6 CS ULTIMA 1/2 Circle   Diam. 150 microns	30cm Black		U7000				
TG160-6 ADVANCED MICRO POINT 1/2 Circle   Diam. 150 microns	30cm Black		W1756				
6.5mm							
BV130-5 3/8 Circle   Diam. 130 microns chord length 5.1 mm	13cm Black		W2810		W2808		
BV130-5 3/8 Circle   Diam. 130 microns chord length 5.1 mm	30cm Black		W2801				
CS140-6 CS ULTIMA 3/8 Circle   Diam. 150 microns	30cm Black		U7003				
TG100-8 MICRO POINT 1/4 Circle   Diam. 200 microns chord length 5.9 mm	30cm Black				1714G		
TG140-6 ADVANCED MICRO POINT 3/8 Circle   Diam. 150 microns	30cm Black	W1780	W1770	W1769			
TG140-8 MICRO POINT 3/8 Circle   Diam. 200 microns	30cm Black		W1718				

NON-ABSORBABLE
SYNTHETIC
MONOFILAMENT







(Polyamid 6 and 66)Suture **ETHILON**

Needle and Description	Length Colour	Size / Product Code						
		10/0	9/0	8/0	7/0	6/0	5/0	4/0
8mm								
G-7 MICRO POINT 1/2 Circle ▽	45cm Black			W1765				
P-6 PRIME 3/8 Circle ▽	45cm Black				1647G	W1600T	W1601T	
9.3mm								
BV-1 3/8 Circle ⊙	45cm Black			W8170				
11mm								
P-1 PRIME 3/8 Circle ▽	45cm Black				1696G	W1610T	W1611T	
	45cm Clear					689G		
	45cm Blue							W1612T
13mm								
P-3 MULTIPASS PRIME 3/8 Circle ▽	45cm Clear						690G	
16mm								
CFS-3 3/8 Circle △	45cm Black					W511H		
PC-3 PRIME 3/8 Circle △	45cm Blue						W1616T W526	
	45cm Black					W507		

ETHILON

ETHILON (Polyamid 6 and 66) Suture















NON-ABSORBABLE
SYNTHETIC
MONOFILAMENT

Needle and Description	Length Colour	Size / Product Code						
		6/0	5/0	4/0	3/0	2/0	0	1
19mm								
PC-5 PRIME 3/8 Circle 	45cm Blue			W519				
FS-2 3/8 Circle 	45cm Blue			W319				
PS-2 PRIME 3/8 Circle 	45cm Blue			W1620T	W1621T			
	45cm Clear			1611G				
26mm								
FS 3/8 Circle 	45cm Blue				W320			
PS PRIME 3/8 Circle 	45cm Blue				W1625T			
	75cm Blue				W1685T	W1626T		
31mm								
MH-1 1/2 Circle 	100cm Black					W568		

ETHILON

NON-ABSORBABLE
SYNTHETIC
MONOFILAMENT





(Polyamid 6 and 66)Suture **ETHILON**

Needle and Description	Length Colour	Size / Product Code						
		6/0	5/0	4/0	3/0	2/0	0	1
40mm								
CT <small>1/2 Circle</small>  	100cm Black						W746	W747
CEC135-40 <small>3/8 Circle</small>  	45cm Blue						W786	
	75cm Blue						W1632T	
CP <small>1/2 Circle</small>  	100cm Black						W737	
45mm								
LS-1 <small>3/8 Circle</small>  	100cm Black						W736	
48mm								
CTX <small>1/2 Circle</small>  	100cm Black							W749
CPX <small>1/2 Circle</small>  	100cm Black							W738
60mm								
KS STRAIGHT  	75cm Blue						W788	

ETHILON

ETHILON (Polyamid 6 and 66) Suture

NON-ABSORBABLE
SYNTHETIC
MONOFILAMENT

Needle and Description	Length Colour	Size / Product Code				
		0	1	2	3	4
90mm						
CKP-3 3/8 Circle  	100cm Blue	W796	W797	W798		
CKP-3 3/8 Circle   <p>Each with 40mm retention tubing</p>	2x50cm Black		W2797			