



PROLENE (Polypropylene) Suture

CLASSIFICATION

Non-Absorbable(비흡수성)

조직을 강하고 안정적으로 지지하면서 감염된 부위에서도 강도가 약해지지 않고 변함없는 인장강도를 유지한다.

Synthetic(합성사)

Polypropylene으로 만들어져 Tissue Reaction이 최소화된다.

Monofilament(단선사)

표면이 매끄러워 조직 투과력이 뛰어나고 박테리아의 증식을 최소화 할 수 있어 감염의 위험이 적다.

PLIABILITY

유연성

단선사이지만 매우 유연하여 실의 꼬임(Memory) 현상이 적고 매듭이 쉬우며 안정적이다.

COLOUR

색상

파란색


















INDICATION & CHARACTER

주로 흉부외과에서 주의를 요하는 조직에 사용 할 수 있으며, 혈전화를 일으키지 않기 때문에 혈관 수술에 있어서 최적의 봉합사이다.

특히 Visi-Black, CC, Hemo-Seal 등 특수한 형태의 봉합침이 달려 있는 제품들이 있다.

NON-ABSORBABLE
SYNTHETIC
MONOFILAMENT

(Polypropylene)Suture **PROLENE**

Needle and Description	Length Colour	Size / Product Code						
		10/0	9/0	8/0	7/0	6/0	5/0	4/0
3.8mm								
BV75-3 <small>38 blade</small> 	 Diam. 75 microns chord length 3.0 mm	13cm Blue	W2790					
4.8mm								
CSC-6 COMPOUND CURVE <small>38 blade</small> 		30cm Blue	7796G					
5mm								
BV130-4 <small>38 blade</small> 	 Diam. 130 microns chord length 3.7 mm	13cm Blue		W2780	W2779			
5.1mm								
BV100-4 <small>38 blade</small> 	 Diam. 100 microns chord length 4.0 mm	13cm Blue	W2794	W2783				
5.5mm								
CS160-6 <small>102 blade</small> 		30cm Blue	7761G					
5.5mm / 13mm								
CS160-6 <small>ULTIMA 1/2 circle</small> 	 Diam. 150 microns	10cm Blue	9090G					
CTC-6 <small>ULTIMA 1/4 circle</small> 								
6.5mm								
BV130-5 <small>38 blade</small> 	 Diam. 130 microns chord length 5.1 mm	13cm Blue			W2775			
BV130-5 MULTIPASS <small>38 blade</small> 	 Diam. 130 microns chord length 5.1 mm	45cm Blue			W2777			
		60cm Blue			8732H			

PROLENE

PROLENE (Polypropylene) Suture

NON-ABSORBABLE
SYNTHETIC
MONOFILAMENT

Needle and Description	Length Colour	Size / Product Code							
		10/0	9/0	8/0	7/0	6/0	5/0	4/0	
6.5mm									
TG140-6 MICROPOINT 3/8 circle 	 Diam. 150 microns	30cm Blue	W1777						
TG140-8 MICROPOINT 3/8 circle 	 Diam. 200 microns	23cm Blue		W1709					
TG140-8 MICROPOINT 3/8 circle 	 Diam. 200 microns	30cm Blue	W1710						
8mm									
BV175-6 MULTIPASS 3/8 Circle 		45cm Blue			8740H	8734H			
		60cm Blue			8741H	8735H			
		2x75cm Blue				M8733			
CC175-6 3/8 Circle 		60cm Blue			W8101 VB				
P-6 PRIME 3/8 Circle 		45cm Blue					W8000T		
9.3mm									
BV-1 3/8 Circle 		45cm Blue				8701H 8301H VB	8806H		
		60cm Blue				8702H W8304 VB	W8712 W8305 VB	9702H HS	
		75cm Blue				8703H	8709H		
		4x75cm Blue				M8703			








PROLENE

VB VISI-BLACK

HS HEMO-SEAL

NON-ABSORBABLE
SYNTHETIC
MONOFILAMENT

(Polypropylene)Suture **PROLENE**






Needle and Description	Length Colour	Size / Product Code						
		10/0	9/0	8/0	7/0	6/0	5/0	4/0
9.3mm								
BV175-8 MULTIPASS 3/8 Circle 	45cm Blue			8752H	8746H			
	60cm Blue			8753H	W8702			
	75cm Blue			8754H	8755H			
CC MULTIPASS 3/8 Circle 	45cm Blue				8705H			
	60cm Blue			W8703	W8704	W8807	W8662 HS	
11mm								
BV 3/8 Circle 	60cm Blue					W8597		
	75cm Blue					R776H		
RB-3 1/2 Circle 	60cm Blue				8206H			
CC-5 1/2 Circle 	60cm Blue						W8803	
CC-11 MULTIPASS 3/8 Circle 	60cm Blue				W8801	W8802		
P-1 PRIME 3/8 Circle 	45cm Blue					W8697 W8003T		

PROLENE

HS HEMO-SEAL

PROLENE (Polypropylene) Suture










NON-ABSORBABLE
SYNTHETIC
MONOFILAMENT

Needle and Description	Length Colour	Size / Product Code						
		8/0	7/0	6/0	5/0	4/0	3/0	2/0
17mm								
RB-1 1/2 Circle 	75cm Blue				8870H	8871H	8872H	
RB-1 1/2 Circle 	45cm Blue				8756H	8757H		
	90cm Blue				W8556	W8557	W8558	
BB 3/8 Circle 	90cm Blue				8580H	8581H		
V-5 TAPERCUT 1/2 Circle 	90cm Blue					W8935	W8936	W8937
19mm								
PC-5 PRIME 3/8 Circle 	45cm Blue					W8391	W8392	

PROLENE

NON-ABSORBABLE
SYNTHETIC
MONOFILAMENT

(Polypropylene)Suture **PROLENE**













Needle and Description	Length Colour	Size / Product Code							
		5/0	4/0	3/0	2/0	0	1	2	
20mm									
SH-2 1/2 Circle ⊙		90cm Blue		W8761 W8340 VB					
CC-20 1/2 Circle ⊙		90cm Blue		W8840					
22mm									
SH-1 1/2 Circle ⊙		75cm Blue			W8770				
SH-1 1/2 Circle ⊙		75cm Blue		8526H					
26mm									
SH 1/2 Circle ⊙		75cm Blue		8831H	8832H	8833H	8834H		
SH 1/2 Circle ⊙		90cm Blue		8521H	W8522	8523H			
		120cm Blue			8534H	8533H			
CC-25 1/2 Circle ⊙		90cm Blue		W8845					
V-7 TAPERCUT 1/2 Circle ⊙		90cm Blue			W8976	W8977			
PC-25 PRIME 3/8 Circle △		45cm Blue			W538				

PROLENE

VB VISI-BLACK

PROLENE (Polypropylene) Suture













NON-ABSORBABLE
SYNTHETIC
MONOFILAMENT

Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	0	1	2
31mm								
MH-1 1/2 Circle  	75cm Blue					W295		
MH-1 1/2 Circle  	90cm Blue			W8525				
MO-5 1/2 Circle  	100cm Blue					W8430		
36mm								
CT-1 1/2 Circle  	75cm Blue			8422H	8423H	8424H	8425H	
MH 1/2 Circle  	90cm Blue			8842H	8843H			
	120cm Blue					8853H		
V-34 TAPERCUT 1/2 Circle  	100cm Blue					W486		

PROLENE

NON-ABSORBABLE
SYNTHETIC
MONOFILAMENT





(Polypropylene)Suture **PROLENE**

Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	0	1	2
40mm								
CT <small>1/2 Circle</small>  	75cm Blue				8433H	8434H		
	100cm Blue						W742	
MO-2 <small>1/2 Circle</small>   <p>Heavier wire diameters</p>	100cm Blue					W8434		
CC-40 <small>1/2 Circle</small>  	90cm Blue			W8851				
CP <small>1/2 Circle</small>  	100cm Blue					W8470		
CFLX <small>3/8 circle</small>  	100cm Blue				W8689			
48mm								
CPX <small>1/2 Circle</small>  	100cm Blue						W8475	

PROLENE

PROLENE (Polypropylene) Suture

NON-ABSORBABLE
SYNTHETIC
MONOFILAMENT

Needle and Description	Length Colour	Size / Product Code				
		3/0	2/0	0	1	2
60mm						
KS STRAIGHT △ 	75cm Blue		W8623			
70mm						
ST-70 STRAIGHT ⊙ 	75cm Blue		W8400			
75mm						
KS-75 STRAIGHT ▽ 	100cm Blue	W751	W752	W753		
90mm						
CKP-3 3/8 Circle △  Each with 40mm retention tubing	2x50cm Blue				W2989	

PROLENE