



ETHIBOND (Polybutylate Coated Polyester) Suture

CLASSIFICATION

Non-Absorbable(비흡수성)

조직내에서 영구적인 인장강도를 유지한다.

Synthetic(합성사)

Polyester 재질로 Tissue Reaction이 최소화된다.

Braided(복선사)

강도가 매우 뛰어나고 핸들링 및 매듭의 안정성이 우수하다.

COATING

코팅재질

원재료(Polyester)와 결합력이 뛰어난 Polybutylate로 코팅되어 있어 코팅재질이 잘 벗겨지지 않아 강도가 아주 뛰어나고 조직 통과시 조직에 손상을 유발하지 않는다.

COLOUR

색상






녹색 또는 흰색

INDICATION & CHARACTER

강도가 뛰어나고 핸들링 및 매듭의 안정성이 우수하여 주로 심장 판막 수술 및 정형외과, 안과 수술에 많이 사용된다. 특히 심장용 봉합사의 경우, 녹색과 흰색 두가지가 있어서 한 번에 많은 봉합사를 사용할 때 구분되어 편리하다.















NON-ABSORBABLE
SYNTHETIC
BRAIDED

(Polybutylate Coated Polyester)Suture **ETHIBOND**

Needle and Description	Length Colour	Size / Product Code						
		7/0	6/0	5/0	4/0	3/0	2/0	0
8mm								
S-14 1/4 Circle 	45cm White			W893				
	45cm Green			W894				
S-22 1/2 Circle 	45cm White			W883				
11mm								
BV 3/8 Circle 	75cm Green		W6597					
13mm								
C-1 3/8 Circle 	75cm Green		W6890	W6891				
17mm								
RB-1 1/2 Circle 	60cm Green		W6757		W6759			
	90cm Green						W6760	







ETHIBOND

ETHIBOND (Polybutylate Coated Polyester) Suture

Needle and Description	Length Colour	Size / Product Code						
		5/0	4/0	3/0	2/0	0	1	2
17mm								
V-5 TAPERCUT 1/2 circle  	90cm Green	W6934	W6935	W6936	W6937			
	90cm White				W6917			
V-5 VALVE KIT TAPERCUT 1/2 circle  	10x75cm 5 White 5 Green				W10B52			
V-5 VALVE KIT TAPERCUT 1/2 circle   <p>Firm pledget 3mm x 3mm x 1.5mm</p>	10x75cm 5 White 5 Green				W10B54			
V-5 VALVE KIT TAPERCUT 1/2 circle   <p>Firm pledget 6mm x 3mm x 1.5mm</p>	10x75cm 5 White 5 Green				W10B55			
20mm								
SH-2 1/2 Circle  	75cm Green		W6761					
22mm								
SH-1 VALVE KIT 1/2 Circle   <p>Firm pledget 7mm x 3mm x 1.5mm</p>	10x75cm 5 White 5 Green				PXX92			
OS-4 1/2 Circle  	75cm Green							X519H









NON-ABSORBABLE
SYNTHETIC
BRAIDED

(Polybutylate Coated Polyester) Suture **ETHIBOND**

Needle and Description	Length Colour	Size / Product Code						
		4/0	3/0	2/0	0	1	2	3
26mm								
SH 1/2 Circle 	75cm Green	W6831	W6832					
SH 1/2 Circle 	90cm Green			W6767				
	100cm Green		W6552					
V-7 TAPERCUT 1/2 Circle 	90cm Green		W6976	W6977	W6978			
	90cm White			W6987				
V-7 VALVE KIT TAPERCUT 1/2 circle 	10x75cm 5 White 5 Green			W10872				
V-7 VALVE KIT TAPERCUT 1/2 circle  Firm pledget 6mm x 3mm x 1.5mm	10x75cm 5 White 5 Green			W10877				
31mm								
MH-1 1/2 Circle 	75cm Green				W975			

ETHIBOND

ETHIBOND (Polybutylate Coated Polyester) Suture

Needle and Description	Length Colour	Size / Product Code				
		1	2	3	4	5
40mm						
OS-8 1/2 Circle  	75cm Green	X548H				
45mm						
V-39 TAPERCUT HEAVY 1/2 Circle  	4x75cm Green		W4843			
55mm						
V-55 TAPERCUT HEAVY 1/2 Circle  	4x75cm Green				W4846	
90mm						
CKP-3 3/8 Circle   <p>Each with 40mm retention tubing</p>	2x50cm Green	W2991				

NON-ABSORBABLE
SYNTHETIC
BRAIDED

ETHIBOND NON-NEEDED

Description	Length Colour	Size / Product Code						
		4/0	3/0	2/0	0	1	2	3
SUTUPAK* Pre-Cut Sutures for Ligation	13x60cm Green		W6232	W6233	W6234			
Sutures for Ligation	180cm Green					W6155	W6156	